

ID	Requirement	Related Use Case(s)	Implemented By	Tested By	Test Description
REQ-01	The system shall allow users to create new profiles with custom settings	Use Case 1: Create a New Profile	<ul style="list-style-type: none"> <li>- Simulation::createProfile()</li> <li>- InsulinPump constructor</li> </ul>	<ul style="list-style-type: none"> <li>- ProfileCreationTest</li> <li>- IntegrationTest_ProfileManagement</li> </ul>	<ul style="list-style-type: none"> <li>- Verify profile is created with specified name</li> <li>- Verify default settings are applied</li> <li>- Test duplicate name handling</li> </ul>
REQ-02	The system shall allow users to delete existing profiles	Use Case 2: Delete a Profile	<ul style="list-style-type: none"> <li>- Simulation::deleteProfile()</li> </ul>	<ul style="list-style-type: none"> <li>- ProfileDeletionTest</li> <li>- IntegrationTest_ProfileManagement</li> </ul>	<ul style="list-style-type: none"> <li>- Verify profile is removed from system</li> <li>- Verify confirmation before deletion</li> <li>- Test handling when current profile is deleted</li> </ul>
REQ-03	The system shall allow users to switch between profiles	Use Case 1: Create a New Profile	<ul style="list-style-type: none"> <li>- Simulation::switchProfile()</li> </ul>	<ul style="list-style-type: none"> <li>- ProfileSwitchingTest</li> <li>- IntegrationTest_ProfileManagement</li> </ul>	<ul style="list-style-type: none"> <li>- Verify active profile changes</li> <li>- Verify settings update to reflect new profile</li> <li>- Test invalid profile selection</li> </ul>
REQ-04	The system shall calculate and deliver manual bolus insulin based on glucose level and carbohydrate intake	Use Case 3: Deliver a Manual Bolus	<ul style="list-style-type: none"> <li>- InsulinPump::manualBolus()</li> <li>- InsulinPump::deliverBolus()</li> </ul>	<ul style="list-style-type: none"> <li>- ManualBolusTest</li> <li>- IntegrationTest_InsulinDelivery</li> </ul>	<ul style="list-style-type: none"> <li>- Verify bolus calculation based on inputs</li> <li>- Verify insulin reservoir is updated</li> <li>- Test delivery with insufficient insulin</li> </ul>
REQ-05	The system shall deliver basal insulin automatically at configured rates	Use Case 4: Automatically Deliver Basal Insulin	<ul style="list-style-type: none"> <li>- InsulinPump::setBasalRates()</li> <li>- Simulation::updateSimulation()</li> </ul>	<ul style="list-style-type: none"> <li>- BasalDeliveryTest</li> <li>- IntegrationTest_InsulinDelivery</li> </ul>	<ul style="list-style-type: none"> <li>- Verify basal insulin delivery at scheduled intervals</li> <li>- Verify correct amounts delivered</li> <li>- Test basal rate changes</li> </ul>

REQ-06	The system shall maintain a historical record of all insulin delivery events	Use Case 5: View Insulin Pump History	<ul style="list-style-type: none"> <li>- InsulinPump::logBolusEvent()</li> <li>- InsulinPump::getHistory()</li> <li>- InsulinPump::printRecords()</li> </ul>	<ul style="list-style-type: none"> <li>- HistoryLoggingTest</li> <li>- IntegrationTest_DataStorage</li> </ul>	<ul style="list-style-type: none"> <li>- Verify events are logged accurately</li> <li>- Verify all required data is captured</li> <li>- Test history retrieval and formatting</li> </ul>
REQ-07	The system shall allow users to view the insulin pump history	Use Case 5: View Insulin Pump History	<ul style="list-style-type: none"> <li>- Simulation::storeCurrentProfileRecords()</li> <li>- Simulation::printAllLogs()</li> </ul>	<ul style="list-style-type: none"> <li>- HistoryViewingTest</li> <li>- IntegrationTest_DataStorage</li> </ul>	<ul style="list-style-type: none"> <li>- Verify history display functionality</li> <li>- Verify filtering and sorting options</li> <li>- Test with large history datasets</li> </ul>
REQ-08	The system shall allow users to suspend and resume insulin delivery	Use Case 6: Suspend/Resume Insulin Delivery	<ul style="list-style-type: none"> <li>- InsulinPump::suspendInsulinDelivery()</li> <li>- InsulinPump::resumeInsulinDelivery()</li> </ul>	<ul style="list-style-type: none"> <li>- SuspendResumeTest</li> <li>- IntegrationTest_SafetyFeatures</li> </ul>	<ul style="list-style-type: none"> <li>- Verify insulin delivery stops when suspended</li> <li>- Verify delivery resumes properly</li> <li>- Test bolus attempts during suspension</li> </ul>
REQ-09	The system shall provide PIN security for critical settings	Use Case 3: Deliver a Manual Bolus Use Case 6: Suspend/Resume Insulin Delivery	<ul style="list-style-type: none"> <li>- InsulinPump::password()</li> <li>- InsulinPump::setPIN()</li> </ul>	<ul style="list-style-type: none"> <li>- PINSecurityTest</li> <li>- IntegrationTest_SafetyFeatures</li> </ul>	<ul style="list-style-type: none"> <li>- Verify PIN validation for sensitive operations</li> <li>- Test incorrect PIN handling</li> <li>- Verify PIN change functionality</li> </ul>
REQ-10	The system shall monitor and report battery status	Use Case 7: Battery Status Check	<ul style="list-style-type: none"> <li>- InsulinPump::checkBatteryStatus()</li> </ul>	<ul style="list-style-type: none"> <li>- BatteryMonitoringTest</li> <li>- IntegrationTest_SafetyFeatures</li> </ul>	<ul style="list-style-type: none"> <li>- Verify low battery alerts</li> <li>- Test critical battery threshold actions</li> <li>- Verify battery level display accuracy</li> </ul>

REQ-11	The system shall check insulin levels and alert when low	Use Case 3: Deliver a Manual Bolus Use Case 4: Automatically Deliver Basal Insulin	<ul style="list-style-type: none"> <li>- InsulinPump::checkInsulinLevel()</li> <li>- InsulinPump delivery methods</li> </ul>	<ul style="list-style-type: none"> <li>- InsulinLevelTest</li> <li>- IntegrationTest_SafetyFeatures</li> </ul>	<ul style="list-style-type: none"> <li>- Verify low insulin alerts</li> <li>- Test delivery attempts with low insulin</li> <li>- Verify proper threshold detection</li> </ul>
REQ-12	The system shall automatically adjust insulin delivery based on CGM data (Control-IQ)	Use Case 4: Automatically Deliver Basal Insulin	<ul style="list-style-type: none"> <li>- InsulinPump::integrateCGM()</li> <li>- Control-IQ algorithm</li> </ul>	<ul style="list-style-type: none"> <li>- ControlIQTest</li> <li>- IntegrationTest_AutomaticAdjustment</li> </ul>	<ul style="list-style-type: none"> <li>- Verify automatic suspension for low glucose</li> <li>- Test increased delivery for high glucose</li> <li>- Verify proper response to glucose trends</li> </ul>
REQ-13	The system shall detect and handle pump malfunctions and errors	All Use Cases	<ul style="list-style-type: none"> <li>- Error handling in all components</li> <li>- InsulinPump::logError()</li> </ul>	<ul style="list-style-type: none"> <li>- ErrorHandlingTest</li> <li>- IntegrationTest_SafetyFeatures</li> </ul>	<ul style="list-style-type: none"> <li>- Verify error detection for various scenarios</li> <li>- Test appropriate alert generation</li> <li>- Verify system recovery after errors</li> </ul>
REQ-14	The system shall provide data visualization for glucose and insulin trends	Use Case 5: View Insulin Pump History	<ul style="list-style-type: none"> <li>- Data visualization components</li> <li>- Reporting functionality</li> </ul>	<ul style="list-style-type: none"> <li>- VisualizationTest</li> <li>- IntegrationTest_Reporting</li> </ul>	<ul style="list-style-type: none"> <li>- Verify graph generation for time-series data</li> <li>- Test visualization of different metrics</li> <li>- Verify interactive data exploration</li> </ul>
REQ-15	The system shall support extended bolus delivery over time	Use Case 3: Deliver a Manual Bolus	<ul style="list-style-type: none"> <li>- InsulinPump::deliverExtendedBolus()</li> </ul>	<ul style="list-style-type: none"> <li>- ExtendedBolusTest</li> <li>- IntegrationTest_InsulinDelivery</li> </ul>	<ul style="list-style-type: none"> <li>- Verify immediate portion delivery</li> <li>- Test scheduled portion delivery</li> <li>- Verify proper scheduling and timing</li> </ul>